

| Test No. | Product          | Weight % | Assays, g/t, % |        |      |       |       | % Distribution |        |
|----------|------------------|----------|----------------|--------|------|-------|-------|----------------|--------|
|          |                  |          | Au             | Ag     | Fe   | Sb    | As    | Au             | Ag     |
| 4        | Cu Cl. Conc.     | 1.35     | 8.32           | 11675  | -    | -     | -     | 11.99          | 26.01  |
|          | Pb Cl. Conc.     | 5.56     | 6.32           | 2691   | 8.61 | 4.82  | 0.66  | 37.60          | 24.75  |
|          | Cu-Pb Ro. Conc.  | 18.61    | 4.01           | 3165.5 | -    | -     | -     | 79.77          | 97.36  |
|          | Zn Cl. Conc.     | 5.00     | 0.49           | 67.8   | 2.80 | 0.026 | 0.032 | 2.62           | 0.56   |
|          | Zn Ro. Conc.     | 8.79     | 0.67           | 75.89  | 3.20 | 0.05  | 0.04  | 6.25           | 1.10   |
|          | Bulk Ro. Conc.   | 27.40    | 2.94           | 2174.5 | -    | -     | -     | 86.03          | 98.46  |
|          | Zn Ro. Tail.     | 72.60    | 0.18           | 12.80  | -    | -     | -     | 13.97          | 1.54   |
|          | Head (Calc.)     | 100.00   | 0.94           | 605.04 | -    | -     | -     | 100.00         | 100.00 |
| 6        | Pb Cl. Conc.     | 7.01     | 5.96           | 5375   | 4.25 | 9.37  | 1.27  | 61.44          | 60.25  |
|          | Pb Ro. Conc.     | 19.96    | 2.64           | 3037.6 | -    | -     | -     | 77.41          | 96.91  |
|          | Zn Cl. Conc.     | 3.18     | 0.94           | 174.0  | 7.13 | 0.10  | 0.07  | 4.40           | 0.88   |
|          | Zn Ro. Conc.     | 15.57    | 0.45           | 78.51  | -    | -     | -     | 10.27          | 1.95   |
|          | Pb-Zn Bulk Conc. | 35.53    | 1.68           | 1740.9 | 6.42 | 2.92  | 0.42  | 87.68          | 98.87  |
|          | Zn Ro. Tail.     | 64.47    | 0.13           | 11.0   | -    | -     | -     | 12.32          | 1.13   |
|          | Head (Calc.)     | 100.00   | 0.68           | 625.60 | -    | -     | -     | 100.00         | 100.00 |
| 7        | Cu Cl. Conc.     | 1.36     | 7.71           | 590.0  | -    | -     | -     | 12.63          | 13.20  |
|          | Pb Cl. Conc.     | 5.56     | 6.06           | 4248.0 | 8.62 | 7.11  | 1.01  | 40.60          | 38.79  |
|          | Cu-Pb Ro. Conc.  | 13.80    | 4.90           | 4070.8 | -    | -     | -     | 81.49          | 92.33  |
|          | Zn Cl. Conc.     | 6.74     | 0.47           | 127.0  | 2.92 | 0.12  | 0.045 | 3.82           | 1.41   |
|          | Zn Ro. Conc.     | 23.87    | 0.30           | 157.03 | -    | -     | -     | 8.74           | 6.16   |
|          | Bulk Ro. Conc.   | 37.67    | 1.99           | 1590.7 | 5.96 | 2.64  | 0.35  | 90.23          | 98.49  |
|          | Zn Ro. Tail.     | 62.33    | 0.13           | 14.7   | -    | -     | -     | 9.77           | 1.51   |
|          | Head (Calc.)     | 100.00   | 0.83           | 608.37 | -    | -     | -     | 100.00         | 100.00 |